



STAR MED



TaeWOONG
NITI-S



STAR MED

MANUFACTURER:

TaeWoong Medical Co., Ltd
14, Gojeong-ro, Wolgot-myeon, Gimpo-si, Gyeonggi-do, South Korea 10022
TEL: +82 31 996 0641~4 FAX: +82 31 996 0645~6


DISTRIBUTOR:

STAR MED America
4801 Wilshire Blvd, Suite 101
Los Angeles, CA 90010, USA
E: info@starmed-america.com T: +1.323.692.0028
www.starmed-america.com




©2026 STAR MED America, Inc. All Rights Reserved.
For distribution in the USA only.
SA-MK-BMM-004D(rev 2)

COMPETITOR



Low Conformability

TAEWOONG



High Conformability

Niti-S

Biliary Stent

PRODUCT FEATURES

► Unfixed Cells with Weaving Construction

Unfixed cells resemble a “chain link fence” allowing the ability for the Niti-S Biliary stent to retain its shape

► High Conformability

Optimal combination of radial and axial force to maintain luminal patency in tortuous anatomy

► Precise Placement and Control


Minimal foreshortening for accurate positioning

► Radiopaque Markers

3 (three) at both ends and 2 (two) in the middle

Percutaneous Approach

Item Code	Stent		Delivery	
	Diameter (mm)	Length (cm)	Profile (Fr)	Usable Length (cm)
TD0804	8	4	8	50
TD0806		6		
TD0808		8		
TD0810		10		
TD1004	10	4		
TD1006		6		
TD1008		8		
TD1010		10		



CLINICAL ARTICLE: A comparison of the Niti-D biliary uncovered stent and the uncovered Wallstent in malignant biliary obstruction by Ki Young Yang, MD et al [Gastrointest Endosc. 2009 Jul;70(1):45-51]

INDICATION: The Niti-S Biliary Stent is indicated for the palliation of malignant strictures in the biliary tree.